



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Autoclave YESON E series 22I (New)**

New,

Modern Yeson E-series autoclave with a capacity of 22 litres,

Pre-vacuum fractionated (class B - medical) allows to sterilize a variety of unpackaged and packaged instruments - solid products, porous products, objects with narrow clearance according to classification PN-EN 13060-2015,

Designed for:

doctor's surgeries,  
doctor's offices, dentist's

cosmetics,  
veterinary surgeries,  
clinics,

hospitals,

Features and benefits of E-series autoclave:

Stainless steel chamber,

Class B fractionated vacuum,

Built-in printer,

USB connector for pen drive,

LCD display,

Touch control buttons,

Polish and English menu,

Anti-accidental opening lock,

Automatic fuses,

9 operating programmes,

2 test programmes,

Delayed start function,

Pre-heating function,

Password protection function

Standard equipment:

Enlarged basket with 3 trays,

Power cord,

Tray holder,

Drain hose,

Polish manual,

English instructions,

Technical data:

Capacities: 22l,

Power: 2.3 kW or 2.4 Kw respectively,

Maximum temperature 137° C,

Unit dimensions: 510mm x 680mm x 430mm,

Dimensions of the chamber (22l)  $\Phi$  250 x 450 mm,

Weight: 54 kg,

European CE Certificate,

It has a valid technical passport,

Warranty: 24 months

