



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Stand with heater SCHMITZ MTC10 242.131 (Reconditioned)

Used, perfect technical condition,
German production,
Power supply: 230V,
Frequency: 50/60 Hz,
Power consumption: 400 W,
Heater among others for saline, etc.,
Temperature regulation in the range: 35 - 60 °,
Heating plate covered with Teflon for easy cleaning,
The construction is made of stainless steel,
Manual height adjustment using a gas spring in the range: 101-131 cm,
Stainless steel bowl 18/10 (NOT INCLUDED),
Mobile with wheels, two with brakes,
Dimensions: 60 x 60 x 101/131 cm,
Has a current Technical Passport issued,
1-month warranty

