



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS CYF-4 cystofiberoscope (Reconditioned)

Used, perfect technical condition,
German production,
The camera has 3 small blind spots,
Direction of view: forward view,
Outer diameter: 5.2 mm,
Working channel diameter: 2.4 mm
Working length: 380 mm,
Overall length: 670 mm,
Depth of field: 3-50 mm,
Max. The angle of bending the tip of the endoscope up and down: 210° / 120°
Maximum angle of field of view: 120°,
Included:
Suitcase,
Has a current review issued, and a Technical Passport,
1-month warranty,
Kat 01

