



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Co2 medical bottle 3 kg. 5 liters (Reconditioned)

Used, very good technical condition,
Polish production,
Capacity: 5 liters,
Steel cylinder,
Cylinder weight: 6 kg,
Filling pressure: 200 bar,
Working pressure: 300 bar,
Has an attestation and valid legalization, min. 8 years,
Dimensions: diameter 13 cm x height 55 cm