



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



---

### **Medical air bottle (AIR) 5 liters (Reconditioned)**

Used, very good technical condition,  
Polish production,  
Capacity: 5 liters,  
Steel cylinder,  
Weight of empty cylinder: about 5 kg,  
Filling pressure: 200 bar,  
Working pressure: 300 bar,  
Has a certificate and valid legalization,  
Overall dimensions: diameter 13 cm, height 55 cm