



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical air bottle (AIR) 10 liters (Reconditioned)

Used, very good technical condition,
Polish production,
Capacity: 10 liters,
Steel cylinder,
Weight of the empty cylinder: about 11 kg,
Filling pressure: 200 bar,
Working pressure: 300 bar,
The bottle is filled with air: 1.6 m³,
Has a certificate and valid legalization,
Dimensions: diameter 13,5 cm, height 90 cm
REVIEW PICTURES BOTTLE

