



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006

kw@praiston.pl



APOLLO rehabilitation table (New)

New,

It has a strengthened steel construction coated with powder paint, standard is a two-part sleep surface,

The table is distinguished by an increased maximum load of 350 kg. Height adjustment of the table is carried out with an electric actuator operated by a manual remote control,

It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available for an extra charge,

Options such as:

Armrests,

Hanger,

Foot remote control (instead of manual remote control),

Technical data:

Length: 200 cm

Width to choose: 65/70/80 cm

Height: 48-98 cm

Headrest angle adjustment: -70° to 45° (by gas spring),

Weight: 120 kg

Maximum load: 350 kg

It has a valid Technical Passport,

12-month warranty.

Possibility of financing (Poland only): instalments, leasing, loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

