



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



### **Diabetic analyzer Care Diagnostica EcoMatic (Reconditioned)**

Used, very good technical condition,  
German production (2013),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 30 VA,  
The analyzer is used to determine the concentration of glucose or lactate in whole blood,  
Measurement method: enzymatic-amperometric,  
Short measurement time and camera performance up to 130 tests per hour ensure quick results,  
Touch, color LCD screen,  
Built-in thermal printer,  
Sample detection: automatic,  
Plate capacity: 50 samples, 1 item 'cito', 1 standard, 3 checks,  
Calibrations: start, cyclical, periodic,  
Built-in connectors: RS232, USB,  
Selection of units: mg / dL or mmol / L,  
Memory: up to 1000 results (with date and time of examination),  
Software: in Polish, English or German,  
The volume of the sample: 20 microliters (whole blood, plasma, serum),  
Sample collection: end-to-end capillary tube + 1.5 / 2.0 ml tube, with reagent,  
Measuring range: Glucose 9 - 900 mg / dl,  
Precision: <3% (at 12 mmol / liter); systematic error <= 2.5%,  
Weight: 6 kg,  
Overall dimensions: 40 x 43 x 24.5 cm,  
Has a current Technical Passport issued,  
1-month warranty

