



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Compressor and Mammal EKOM 2V / M (with box, with membrane dryer) (New)

- New,
- Slovak production,
- Compressor DK50 2V + Mammal,
- **Technical data:**
 - Frequency: 230V / 50 (60) Hz
 - Number of dental sets: 1
 - Control voltage: 24 / 0.1 V / A DC / AC
 - Compressor capacity (working pressure 5-7bar):
 - 108 l.min- 1/5 bar
 - Tank capacity: 25l.
 - Operating mode: S1-100%
 - For power supply: to one dental unit,
 - Oil-free,
 - Volume:
 - ≤ 55 dB [A]
 - Dimensions (W x D x H): 560x684x1251 mm
 - Weight: 126 kg
- Has a valid Technical Passport issued,
- 12 months warranty,
- Financing option (Poland only): Installments, Leases, Loans

