



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SUCTION PUMP FOR EKOM 2.1 (stationary mammal in a box) (New)

- New
- Slovak production,
- Mammal in silent housing size DK50 S,
- **Technical data:**
 - Rated voltage / frequency:
 - 230 V / 50 (60) Hz
 - 110V / 60 Hz
 - Maximum current: 3.3 (3.7) A / 7.2 A
 - Maximum flow: 1100 lit.min - 1
 - Maximum negative pressure: 12 kPa
 - Control voltage: 24V / 0.1A DC / AC
 - Operating mode: S1-100%
 - Volume without box: ≤ 58 dB [A]
 - Volume with box: ≤ 46 dB [A]
 - Weight without box: 15 kg
 - Weight with box: 27 kg
 - Dimensions without box (W x D x H): 290x430x350
 - Dimensions with box (W x D x H) * 380x525x575
- Has a valid Technical Passport issued,
- 12 months warranty,
- Financing option (Poland only): Installments, Leases, Loans

