



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SUCTION PUMP FOR EKOM 2.1Z (stationary mammal without housing) (New)

New
Slovak production,
Mammal without silencing housing,
Mammal comes standard with a control module that turns on the mammal after 24V AC / DC signal flow from the unit,
Technical data:
Rated voltage / frequency:
230 V / 50 (60) Hz
110V / 60 Hz
Max. Current: 3.3 (3.7) A / 7.2 A
Maximum flow: 1100 liters min⁻¹
Maximum vacuum: 12 kPa
Control voltage: 24V / 0.1A DC / AC
Operating mode: S1-100%
Volume without the box: ≤58 dB [A]
Volume with a box: ≤46 dB [A]
Weight without box: 15 kg
Weight with case: 27 kg
Dimensions without box (W x H x W): 290x430x350
Dimensions with the box (sz x g x w) * 380x525x575
Has a current Technical Passport issued,
12 months warranty,
Possibility of financing (only Poland): Installments, Leasing, Loan

