



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Recorder CD-DVD MEDISCO DIVICOM CD VTR-160 (Ex-Display)**

- New, exhibition,
- German production,
- power supply: 230 V,
- Frequency: 50/60 Hz,
- power: 50 W,
- Supported formats:
  - Playback: CD-DVD, HQ-VCD, S-VCD, VCD, audio CD,
  - Recording: CD-DVD, HQ-VCD, S-VCD, VCD,
- Video Format: CCIR 601 PAL
- CD-DVD: 704 x 576 at 25 bps 5 MB/s, 18/20 min,
- S-VCD: 352 x 576 at 25 bps 2.4 MB/s, 36/40 min,
- HQ-VCD: 352 x 288 at 25 bps 2.4 MB/s, 36/40 min,
- VCD: 352 x 288 at 25 bps 1.15 MB/s, 74/80 min,
- In-put Video: S-Video (Y/C mini DIN, Video BNC, DV (IEEE 1394/iLINK) 4-pin connection (CD-VTR 160 only),
- Out-put Video: S-Video (Y/C) mini DIN, Video BNC, Y.U.V., R.G.B., BNC,
- In-put Audio: Stereo cinch,
- Out-put Audio: Stereo cinch,
- **Includes:**
  - DV cord,
  - S-Video,
  - Power cord,
- Dimensions: 220 x 415 x 115 mm,
- Weight: 7 kg,
- User manual in English and German,
- It has exoised the current Technical Passport,
- Warranty:
  - 6 months for national field - Poland,
  - 1 month for international field,
- The possibility of financing (Polish only): Installments, lease, loan.

