



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Recorder DVD SONY DVO1000MD (Reconditioned)

- Used, very good technical condition,
- Japanese production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Recording format:
 - Video: Compression MPEG-2,
 - Audio: Dolby Digital Format,
- Recording medium:
 - DVD+RW
- Recording speed: 2,4x, 4x,
- Thanks to the built-in 80 GB hard disk, recording starts when you press the REC button. What's more, the hard disk also provides additional security for recordings, preventing data loss in the event of an unexpected power outage or damage to the DVD,
- **In-put / Out-put:**
 - Analog composite signal input: BNC x 2, with loop-through,
 - S-Video input: 4-pin DIN x 2, with loop-through,
 - Analog composite output: BNC x 1,
 - S-Video output: 4-pin DIN x 1,
 - Analog audio input: RCA x 2 (L, P),
 - Analog audio output: RCA x 2 (L, P),
 - Monitor audio output: RCA x 1,
 - i.LINK (DV IN): 6-pin i.LINK x 1,
 - RS-232C x 1, 9-pin D-sub,
 - Pilot 1 x 1, mini-jack stereo (for SVRM-100A i DSRM-10 pilots),
 - Pilot 2 x 1, mini jack stereo (for foot switch FS-20),
- Weight: 6 kg,
- Dimensions 212 x 128.5 x 382 mm,
- User manual in English (PDF),
- Cat 02,
- It has exoised the current Technical Passport,
- Warranty:
 - 6 months for national field - Poland,
 - 1 month for international field,
- The possibility of financing (Polish only): Installments, lease, loan.

