



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BTL exercise testing kit (Reconditioned)

Used, very good condition,
UK production,
The set includes:
BTL-08 MT Plus ECG machine;
Used, very good technical condition,
UK production (12.2011),
Power supply: 230 V,
Frequency: 50/60 Hz,
3, 6, 12-lead ECG for printing and recording of 12 leads,
Colour touch screen diag: 5.7" (118 x 89 mm),
Screen resolution: 320 x 240,
Alphanumeric keyboard, function and touch screen buttons,
Contact control of each electrode,
Equipped with thermal printer,
Paper width: 112 mm,
Direct printing from a normal printer (via USB port, independent of PC, printing on A4 office paper),
Sensitivity: 2.5, 5, 10, 20 mm/mV,
Print speed: 5, 10, 25, 50 mm/s,
Line noise filter: 50/60 Hz,
Muscle noise filter: 25, 35 Hz,
Isoline filter: 0.05 (3.2s), 0.11 (1.5s), 0.25 (0.6s), 0.50 (0.3s), 1.50 (0.1s), splines,
Automatic filters: autoadaptive,
Number of printed leads: 3, 4, 6, 12, 3x4+1, 3x4+2, 3x4+3, 4x3+1, 6x2+1, 6x2+2,
Up to 120 examinations in camera memory (10 sec.),
Unlimited number of profiles,
Manual, automatic and LONG mode,
Real-time and synchro recording in automatic mode,
Reverse ECG recording,
Defibrillation protection,
Cardiac pacemaker detection,
Built-in rechargeable battery,
Technical parameters:
Reference frequency: 0.05 - 150 Hz,
Digital resolution: 3.9 µV,
Analogue to digital A/D conversion: 13 bits,
Sampling frequency: ECG 2000 Hz,
Dynamic range: 15.9 mV,
Voltage polarity: ± 400 mV,
Max. voltage compliance: ± 5 V,
Input impedance: > 20 MOhm,
Battery capacity: 3 - 4 hours continuous recording, 15 minutes continuous printing,
Charging time: 10 hours maximum (in case of complete discharge),
Included:
Two 10-lead ECG cables,
Limb clip electrodes: 4 pieces,
Thoracic electrodes: 6 pieces,
Connection cable: RS232 - USB conversion,
USB extension cable,
Power cord,
Instruction manual in Polish,
Dimensions: 33 x 26.5 x 8 cm,
Weight: 3.2 kg,
BTL Treadmill model;
Used, very good technical condition,
UK production (2014),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1.5 KVA,
Computer controlled via RS232 port,
Durable construction with handrails,
Two instantaneous switch-off buttons,
Belt alignment system,
Low step height,
Very quiet function,
Shock-absorbing system,
Special 3 mm thick running belt,
Speed range: 0 - 20 km/h,
Speed step: 0.1 km/h,
Inclination range: 0 - 25 %,
Inclination step: 0.5 %,
Maximum patient weight: 200 kg,
Included:
BTL-08 E300 M exercise test key,
Connection cable,
Power cord,
Effective area: 50 x 150 cm,
Dimensions: 205 x 80 x 112 cm,
Weight: 163 kg,
Comes with a current Technical Passport issued,
Warranty:
6 months for the domestic market - Poland,
1 month for the international market,
Financing options (Poland only): Installments, Leasing, Loan

