



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SURTRON FLASH 160 HF surgical diathermy (New)

New

Italian production,

The SURTRON FLASH 160W HF electrosurgical device is a versatile electrosurgical device that meets the most stringent requirements of the operating theater.

It guarantees patient safety thanks to the possibility of installing a neutral electrode and its constant monitoring.

Output circuit parameters are constantly monitored.

A wide selection of electrodes.

Italian production apparatus.

Possibility of power control ..

Allows constant monitoring of the patient's neutral electrode,

Has functions:

Checking patient circuit

Power control,

Self-tests,

Programs for diathermy:

CUT smooth cut,

ENCHACED coagulation cut,

BLEND mixed cutting,

FORCED COAG forced coagulation,

SOFT COAG soft (white) coagulation without carbonization,

BIPOLAR COAG bipolar coagulation,

BIPOLAR ABL bipolar ablation

Application:

Outpatient, pediatric, vascular, plastic surgery,

Endoscopy,

Dermatology,

Otolaryngology,

Ophthalmology,

Urology,

Dentistry,

Gynecology,

Veterinary medicine,

First aid,

Gastroenterology,

Pneumologia,

Neurosurgery

Equipment:

1 pc bipolar handle with electrode switches (for sterilization),

1 piece. Needle holder for microsurgery,

1 pc. Neutral steel electrode,

Waterproof foot switch on / off mode,

10 pcs of standard electrodes,

5 needles for microsurgery,

3 pcs of 7 cm flat electrodes,

3 needle electrodes 7cm,

3 pcs of 6 cm ball electrodes,

Neutral electrode cable,

Steel neutral electrode,

2 pcs. Disposable neutral electrode,

2 pcs single-sided neutral electrode,

5 m power cable,

Technical data:

Width - 37 cm,

Height - 14.4 cm,

Depth - 28.8 cm,

Curb weight - 6 kg,

Has a valid Technical Passport issued,

12 months warranty

