



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SURTRON ABC surgical diathermy (New)

New

Italian production,

SURTRON is a module that can be connected to a compatible electrosurgical device, specific for electrosurgical coagulation procedures using argon gas, indicated for both open surgery and laparoscopy methods,

Argon is an inert gas which, when used in conjunction with specific high-frequency coagulation currents generated by the electrosurgical unit, allows superficial coagulation without contact due to ionization of the gas flow by high-frequency current,

The SURTRON attachment works in such a way that the connection with a diathermy with compatible properties with argon procedures enables feeding through both manual activation and foot activation. Connection with the SURTRON ABC adapter does not exclude the use of other diathermy functions, in fact, when using compatible LED devices, argon gas supply occurs only when high frequency coagulation (SPRAY)

power is supplied,

Technical data:

Input Voltage: 90-240 V,

Mains frequency: 50/60 Hz,

Electric input power: 70 VA,

Flow regulation: 0.2 - 10 LPM,

Maximum inlet pressure: 3.0 bar (43.5psi),

Electrical class: I CF,

Size WxHxD mm: 470x150x400,

Weight: 8 kg

Has a valid Technical Passport issued,

12 months warranty

