



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CUTERA TITAN lifting laser (Reconditioned)

Used, in very good condition,
of the U.S.A. Production Department,

Power supply: 230 V,

Frequency: 50/60 Hz,

The Cutera Titan laser allows you to non-surgically increase the skin's firmness and achieve a lifting effect,

It is an innovative method of tightening the skin, improving contours of the face, neck, abdomen, thighs and arms. Titan uses far infrared radiation at the boundary of high-power radio waves, it heats the skin to a depth of up to 3 mm,

As a result of radiation, there is a thermal effect in the lower layers of the skin, which causes the collagen fibres to heat up. The heated fibres are irritated and thus stimulated to rebuild. The skin exposed to radiation thickens and is tightened, which results in an improvement of its structure and smoothing of wrinkles,

Safe, comfortable, non-invasive,

For loose and superfluous skin of the whole body,

Stimulation of long-lasting renewal of collagen,

Continuous cooling of the skin during the treatment,

No wear parts (spare parts) - low operating costs,

Wave length: 1100, 1800 nm,

Spot size or treatment area (aiming beam): 10 - 15 mm,

The head of the device is exchanged every 15000 pulses, the current head will allow to perform 9862 more pulses,

He has issued a valid Technical Passport,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (only Poland): Instalments, Leasing, Loan

