



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Reductor for oxygen cylinders with AGA quick coupler (New)

A new one,  
Manufacturer: KORGIEL,  
Model: W8,  
Polish production,  
The regulator is used to lower the pressure of gas taken from the cylinder to the desired value,  
The range of indication of the manometer: 0 - 31,5 MPa,  
Gas inlet pressure: 20 MPa,  
Pressure set on the safety valve: 0.7 MPa,  
Reduced gas pressure:  $0.5 \pm 10\%$  MPa,  
Input connector type: O2 - G3/4",  
Output connector type: AGA quick connector,  
Weight: 0.7 kg,  
The manual in Polish,  
He has a current Technical Passport issued,  
Guarantee until 18.03.2021,  
Possibility of financing (Poland only): Installments, Leasing, Loan

