



## **RESPIRONICS REMSTAR NOVA sleep apnea breathing apparatus (Reconditioned)**

Used, very good condition,

German production,

Power supply: 230 V,

Frequency: 50/60 Hz,

The device with automatic control, which is set to the pressure range and then searches for and provides the lowest pressure level necessary to maintain the airway. Thanks to C-Flex technology, not only does the pressure level change according to needs, but it decreases with every exhalation,

REMstar nova is designed for patients whose pressure needs change at night depending on their sleep stage or position,

Automatic pressure adjustment based on airflow limitations from snoring or apnea events,

The C-Flex technology offers a more comfortable way to carry out sleep therapy by reducing the pressure during exhalation. The degree of pressure reduction depends on one of three settings that you can adjust to your own comfort,

Space for a dedicated memory card for session data recording - no card included,

Works in Auto or CPAP mode,

Automatic starting and stopping of work,

Very quiet operation,

Power cord included,

All placed in the original carrying bag,

The dimensions of the camera: 17.5 x 25 x 12 cm,

The weight of the set: 2.2 kg,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

