



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RESMED VPAP III sleep apnea breathing apparatus (Reconditioned)

Used, very good condition,

Australian production,

Power supply: 230 V,

Frequency: 50/60 Hz,

The device is designed for the treatment of obstructive sleep apnea (OSA) in adult patients in hospital or at home,

Effective VPAP synchronization - patient-machine synchronization is the key to successful treatment. The unique features of Vsync and TiControl VPAP III ensure that the patient is synchronized with the machine even in the event of large leaks and specific respiratory disorders,

Intelligent management of the therapy:

VPAP III offers comprehensive data management functions, enabling your doctor to monitor and manage your therapy,

SmartStart automatically starts and stops the device,

Smart Data provides information on how to use and fit the mask,

The mask fit function allows you to adjust the mask to ensure optimal sealing,

Two modes of operation:

Spontaneous mode,

CPAP mode,

Clear backlit LCD display,

Possibility of connecting a dedicated humidifier - not included,

Very quiet operation,

Comfortable carrying handle,

Power cord included,

All placed in the original carrying bag,

The dimensions of the camera: 23 x 29 x 14.5 cm,

The weight of the set: 3.7 kg,

User manual in English (PDF),

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

