



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Holter PROGETTI PG 8 PLUS 3 channel (New)**

A new one,  
Italian production,  
Designed for the Windows 8® platform, PG-8 Plus includes an exclusive algorithm for impact detection and arrhythmia analysis combined with state-of-the-art digital technology for unparalleled accuracy and reduced scanning time.  
Multiple scanning options: Retrospective, Prospective and Overlay Scanning Modes supported in traditional and advanced scanning options.  
Features and benefits:  
PG-8 Plus software offers all the tools you need to customize your scans for best results,  
Multi View displays for heart rate, ST and ventricular charts combined with full disclosure and full scale ECG,  
A true 3-channel ECG rhythm detection and analysis algorithm provides unparalleled accuracy,  
A real ST analysis,  
Comprehensive morphology (template) editing to add, delete, mark or declare artifacts at the level of individual QRS assemblies or edit an entire morphology class,  
Histograms of events and rhythms,  
Automatic detection of atrial fibrillation,  
Atrial fibrillation and atrial flutter edition,  
Roundabout and 3D charts,  
HRV in time and frequency domain,  
QT analysis,  
Powered by one AAA battery,  
Replaceable SD card,  
128 x 64 dot matrix LCD display,  
Pacemaker analysis,  
Full recording time 24/48 hours,  
Data transfer via Compact Flash card or USB cable,  
Windows 8® / Vista® compatible software,  
Technical data:  
Report formats include  
report summary, heart rate chart, chamber chart and table, ST charts and tables, hourly summaries,  
Automatically detected or selected by the ECG operator in full scale in a variety of bar configurations per page,  
The full test record selected by the site operator,  
Automatic report generation,  
The final report can be completely customized, including in the local language,  
Recommended minimum PC System:  
Intel Pentium IV- 1.3GHz,  
512MB frame (1 GB recommended),  
HD 20 Gbyte,  
Monitor 1024 x 768 graphics resolution,  
Operating System: Windows 8 ®,  
Printer: laser printer,  
Suggested configuration of the PC system:  
Intel Core 2 duo 2.4 GHz,  
RAM 2 GB or more,  
HD 80 Gbyte or more,  
Monitor the SVGA 1024x768 graphic resolution,  
Operating System: Windows 8 ®,  
Printer: laser printer,  
Included:  
3-channel ECG recorder  
USB cable and SD flash card  
Patient cable  
Plastic handle  
Operating instructions  
He has a current Technical Passport issued,  
Warranty: 12 months,  
Possibility of financing: Installments, Leasing, Loan.

