



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **PURITAN BENNETT GOODKNIGHT 420E apnea breathing apparatus (Reconditioned)**

Used, good technical condition,  
French production,  
The power supply will support a 13-volt power supply,  
Frequency: 50/60 Hz,  
The GoodKnight 420 is an auto-adjustable CPAP device that is set to a pressure range and then searches for and provides the lowest pressure level needed to keep the airways open,  
The pressure level will vary as needed, but will not decrease with every exhalation,  
CPAPs with automatic adjustment are designed for patients who cannot participate in the sleep study or whose blood pressure needs change at night depending on their sleep stage or position,  
The GoodKnight 420E model records up to 200 use data sessions (date/time),  
Small in size,  
Simple and intuitive operation,  
Wide pressure range: 4 - 20 cm H<sub>2</sub>O,  
Height compensation function,  
Automatic leakage compensation,  
Adjustable ramp starting pressure,  
The device measures the pressure in the mask, not in the device, for precise measurement and greater comfort,  
It detects and reacts to apnea and snoring,  
Detects and reacts to a restriction in the inspiratory flow,  
It detects apnea with and without heart rate oscillations,  
Very quiet operation,  
Power supply included,  
All placed in the original carrying bag,  
Dimensions: 14.5 x 19.5 x 8 cm,  
The weight of the set: 1.4 kg,  
User manual in English (PDF),  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan

