



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Respiratory assist PAP-device WEINMANN SOMNOSMART 2 (Reconditioned)

- Used,
- Technical condition: very good,
- Visual condition: good,
- Manufactured in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **The SOMNOsmart 2** is a modern apnoea therapy device operating in autoPAP mode. It gently raises the air pressure, preventing apnoea from occurring. By analysing respiratory events such as apnoea or shortness of breath, flow restrictions and snoring, the device automatically adapts the pressure to the needs of the sleeping patient,
- The accuracy of event detection is based on a combination of OPS (oscillation), flow and snoring signals. A unique algorithm guarantees the effectiveness of the therapy,
- The device identifies and differentiates between central and obstructive apnoea,
- The SOMNOsmart 2 is very quiet, easy and comfortable to use (thanks to the AUTO start function and softPAP - pressure reduction in the last phase of exhalation),
- As a sleep apnoea diagnostic device, in combination with the SOMNOsupport software, it provides extensive information presented graphically: RDI, respiratory events, snoring, OPS (Oscillatory Pressure Signal), leaks, etc. and systematises measurement data (therapy calendar),
- Leakage compensation,
- Clear readable display,
- Simple and intuitive operation,
- Very quiet operation: 26 dB,
- **The set includes:**
 - Original carry bag,
 - User manual in multiple languages (PDF),
- Dimensions of the device: 32 x 18.5 x 9 cm,
- Weight of set: 4.6 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Instalments, Leasing, Loan,



