



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



B. BRAUN PERFUSOR COMPACT Infusion (syringe) pump (Reconditioned)

- Used, technical condition very good,
- Manufactured in Germany,
- Power supply: 230 V,
- Emergency power supply: through four new accumulators / batteries LR6 / AA,
- Frequency: 50/60 Hz,
- Power: 10 VA,
- The Perfusor Compact pump can operate for about 10 hours at a flow rate of 99.9 ml/h or about 80 hours at an average flow rate of 10 ml/h,
- Choice of syringes: 20 or 50 ml,
- Delivery rate: from 0.1 to 99.9 ml/h (in 0.1 ml increments),
- Bolus rate: 800 ml/h,
- Infusion volume: from 0.1 to 999.9 ml (in 0.1 ml increments),
- Automatic recognition of syringe size,
- Features foiled control panel and liquid crystal display,
- Audible and visual alarms,
- Multifunctional MFC connection,
- Convenient carrying handle,
- Foldable kit,
- Dimensions: 24.5 x 14.5 x 10.5 cm,
- Weight: 1.4 kg,
- User manual in Polish and English (PDF),
- Cat 04,
- Valid technical passport issued,
- Warranty:
 - 6 months for the domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

