



## **RESPIRONICS BIPAP FOCUS non-invasive ventilation unit (Reconditioned)**

Used, very good condition,  
German production,  
Power supply via 18 V power supply and built-in battery,  
Frequency: 50/60 Hz,  
Respironics BiPAP Focus is a non-invasive hospital grade ventilator designed to treat adult patients with respiratory failure and obstructive apnea during sleep,  
Supports breathing in adult patients weighing 30 kg or more,  
BiPAP Focus offers two ventilation modes:  
Continuous positive output pressure (CPAP), which provides a single positive pressure level for the patient,  
Spontaneous / Temporary (S / T), which provides two levels of positive pressure (one during inhalation and one during exhalation), and performs breaths on time if the patient does not initiate breathing,  
BiPAP system alarms indicate high or low pressure, apnea, patient disconnection, low power or loss of mains power,  
Auto-Trak's digital sensitivity eliminates the need for a perfect seal of the patient interface and constant adjustment, while promoting patient to ventilator synchronization,  
The pressure bar graph shows patient breathing and pressure ranges,  
Simple and intuitive operation,  
Very quiet operation,  
The camera is placed on the original cart:  
Mobile, on four wheels, including two with brakes,  
Spacious basket for accessories and wiring,  
Convenient handle for moving,  
Power supply unit with power cord included,  
The dimensions of the camera: 30.5 x 39 x 21.5 cm,  
Overall dimensions: 50 x 62 x 118 cm,  
The weight of the camera: 4.5 kg,  
The weight of the set: 23 kg,  
The user manual in English,  
He has a current Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
1 month for the international market,  
Possibility of financing (Poland only): Instalments, Leasing, Loan

