



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDELA fat transfer kit (Fat Collector / Fat Transfer) (New)

New,
Swiss production,
The set includes:
1000 ml fat container with a lid and overflow valve,
Reusable stand made of stainless steel,
Reusable silicone hose,
Additional equipment for an extra charge:
Disposable 17 cm hose with clamp and Lauer Lock faucet,
3000 ml fat container
The container is made of high quality, unbreakable, impact-resistant polysulphone. It has a clear marking of liquid capacity. It contains no latex and is suitable for autoclave sterilization,
The cover is made of reinforced polyamide (Santoprene gasket) with a conical patient connector and mechanical overflow protection. It does not contain latex and is suitable for autoclave sterilization,
He has a current Technical Passport issued,
Warranty 12 months,
Possibility of financing (Poland only): Installments, Leasing, Loan.

