



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Docking station for 6 pumps FRESNIUS KABI ORCHESTRA BASE A (Reconditioned)

Used, very good condition,

French production,

Power supply: 230 V or via the built-in 6 V / 1.3 Ah battery,

Frequency: 50/60 Hz,

Power: 125 VA,

The Base A docking station is a data hub and power supply for up to 6 Visio DPS modules and/or Visio DPS modules and/or MVP modules, it allows for faster therapeutic control,

Communication interface compatible with the main PDMS software,

Built-in RS232 interface,

One single power cable for all connected modules,

Battery life - minimum 2 hours,

Equipped with a tripod with a comfortable carrying handle,

No pumps included - it is possible to retrofit the FRESNIUS ORCHESTRA MODULE DPS+ infusion (syringe) pumps for an additional fee,

Dimensions of the docking station:

Height: 55 mm,

Width: 315 mm,

Depth: 200 mm,

Height with tripod: 1015 mm,

The weight of the station with the tripod: 4 kg,

The weight of the docking station with pumps: 24 kg,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

