



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Vital signs patient monitor DATASCOPE DUO (Reconditioned)

Used, very good condition,

U.S.A. production,

Power supply: 230 V and through a 7.4 V / 7800 mAh lithium-ion battery - no battery included (can be retrofitted at extra charge),

Frequency: 50/60 Hz,

The Datascope Duo patient monitor is used for non-invasive and automatic measurement of systolic and diastolic blood pressure, heart rate and oxygen saturation (SpO2) in adults and children,

All blood pressure, heart rate and SpO2 values are shown on clearly legible digital displays,

The device can be used in many different healthcare facilities. These include hospital wards as well as other types of medical facilities such as doctor's offices, clinics and long-term care facilities,

Simple and intuitive operation,

All functions are accessible via the four buttons,

Comfortable carrying handle,

Included:

Air-wire cuff,

Silicone sensor SpO2,

Power cord,

User manual in English (PDF),

All placed on the original tripod:

Basket for wiring and accessories,

Unsymmetrical five-arm base on wheels,

Dimensions of the device:

Width: 13 cm,

Depth: 15 cm,

Height: 24.5 cm,

General dimensions with tripod:

Width: 54 cm,

Depth: 49 cm,

Height: 131.5 cm,

The weight of the set: 12.3 kg,

He has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

