



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Bactericidal lamp UV-C STERILON 108W (New)

New,
Polish production,
UV-C STERILON 108W lamp is designed for surface disinfection,
Using ultraviolet light, it is one of the most effective devices able to remove viruses, bacteria, fungi or destroy the DNA or RNA of any microorganisms that will be exposed,
Simply point the UV lamp directly at the surface for a certain period of time to effectively kill the germs. They are so strong that they also have a harmful effect on human skin or eyes,
Therefore, when using UVC lamps for disinfection, basic safety measures must be observed, i.e. not within the range of the light flux,
The 108W lamps will disinfect the space of about 35m²,
Recommended for use in healthcare facilities, grocery stores, public and private transport, gas stations, courier companies, public administration and education units, as well as for private use in one's own home or car,
Depending on the distance from the light source, temperature and air humidity, we take 15 minutes as an indicative disinfection time,
Features and advantages of UV-C STERILON 108W:
Type of light: UV-C ultraviolet,
Case color: white,
Disinfection efficiency: 35m²,
Lifetime: 8000h,
Power supply voltage: 220-240V,
Degree of tightness: IP20,
Power supply cable 5m,
A multi-positioned arm,
Additional option: ON-OFF switch,
Adapted for mounting on a Future Stand.
Dimensions (A/B/C) [mm]: 270/480/100,
Included:
Lamp,
Universal clip,
Additional accessories (surcharge):
Future Stand,
DEADLINE FOR IMPLEMENTATION OF ORDERS : 48 -72 h (working days) COUNTED FROM THE PAYMENT OF THE PROFORMA INVOICE,
Warranty 12 months

