



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



AGFA CURIX ID identification camera (Reconditioned)

Used, very good condition,

U.S.A. production,

Power supply 230 V,

Frequency 50-60 Hz,

Power: 20 W,

Identification time: 2 s,

Exposure level: 10,

The main feature of the device is the ability to immediately save the date, time and patient data on the membrane in full daylight,

The device is very easy to use,

The total size of the exposed field: 71 x 14 mm,

Size: 39 x 33 x 34 cm.

It has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

1 month for the international market,

Possibility of financing (Poland only): Instalments, Leasing, Loan

