



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006

kw@praiston.pl



Rehabilitation table EOS M02 (New)

New,
 Made in Poland,
 Sturdy steel structure coated with powder paint,
 Standard features include three-section sleep surface, foldable to the Pivot position, made of plywood,
 Height adjustment of the table is carried out with a hand crank at the back of the table frame,
 Fully ball-bearing mounted height adjustment mechanism,
 Individual parts of the tabletop are adjusted by means of gas springs,
 There is a free option to order upholstery in the colour of your choice, except for "PRESTIGE" colours which are available at an extra charge,
 Table features and benefits

Gas spring from the German company BANSBACH,

Upholstery of very good quality,

Foam, 5 cm thick,

Steel ball bearings are used in the movable parts of the table, which makes the table very reliable,

Frame elements made of the best Polish steel,

As additional options are available (for an extra charge):

Armrests,

Side shelves,

Hanger for the underlay,

Castor system,

Euro rails (40cm),

Knee pads,

Injection support,

Technical parameters:

Length 200cm,

Width 65/70/80 cm,

Height 48 - 100 cm,

Headrest angle adjustment -700 to 400,

Adjustment of table top angle 00 to 300,

Own weight 88 kg,

Maximum load 160 kg,

It has a current Technical Passport,

Warranty on construction 24 months,

Warranty on hydraulic accessories 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

