



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Endoscopic table FANOS 01 (New)

New,

Made in Poland,

It has a stable steel construction covered with powder paint,

The table has a column lifting mechanism which eliminates the shift of the tabletop during height adjustment,

In standard, the table has a two-part plywood sleep surface,

Table height adjustment is carried out with an electric actuator operated by a hand-held remote control,

Headrest angle is adjusted manually with a gas spring,

It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available for an extra charge,

Additional options are available (at extra charge):

Armrests,

Paper hanger,

Foot remote control (instead of manual remote control),

Technical data:

Length: 200 cm,

Width: 65/70/80 cm,

Height: 50-90 cm,

Headrest angle adjustment: - 70° to 40°,

Own weight: 88 kg,

Maximum load: 220 kg,

It has a valid Technical Passport,

Warranty on structure 24 months,

Warranty on hydraulic and electrical accessories 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

