



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CAMERON SINO CS-PCL200MX battery (New)

New,
Production P.R.C.,
Battery for the CodeManagment module used to measure CO2 in the Lifepak 20 and Lifepak 20e defibrillators,
Replaces Physio-Control batteries: 11141-000162 / B11827,
Chemical system: Li-Ion,
Rated voltage: 11.1V,
Capacity: 3400 mAh / 37.74 Wh,
Dimensions: 66 x 55 x 21 mm,
Weight: 0.19 kg,
Has a valid Technical Passport issued,
Warranty: 12 months,
Possibility of financing (only Poland): Installments, Leasing, Loan,

