



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SpO2 Adapter przedłużający - GE-Datex Ohmeda (Produkt nowy)

- Nowy,
- Produkcji Urugwajskiej,
- Wtyk do czujnika / wtyk do urządzenia:
 - DB9 / REDEL 10 PIN (BIOOX.DA.IC2),
 - / REDEL 7 PIN (BIOOX.GEO.IC2),
 - DB9 / REDEL 10 PIN (BIOOX.DT.IC2),
 - DB9 / DB9 (BIOOX.GET.IC2),
 - DB9 / BLUE FLAT 11 PIN (BIOOX.GET.IC2),
 - DB9 / HYPERTRONICS 7 PIN (BIOOX.OHT.IC2),
- **Długość kabla 2,0 m** (za dopłatą możliwość zamówienia kabla o dowolnej długości),
- Adapter ma zastosowanie w urządzeniach serii:
 - AS/3, AS/5, CS/3 Modular Monitors with M-series modules, CH/S, CH/2S, CH/RS, CG/2S, Cardiocap I, Cardiocap II, Cardiocap 5, Cardiocap S/5 Light Monitor, Capnomac Ultima, M-Estpr, Oscar II, Satellite, Satellite Plus II, Trans Plus, ULTSVL, ULTSVO, ULTSIO, ULTS00, OSP-200,
 - Ohmeda 3700, 3710, 3770, 3775, 3800, 3900, Biox 3740, 3760, 3800, 4700 Oxicap, 5250 RGM, Sac II, Oxi-Tip Plus, Tuffsat, S/5,
 - CARESCAPE™ modular monitors and S/5, AS/3 and CS/3 modular monitors with M-NE12ST(P)R, M- (N)EST(P)(R) or M-(P)RESTN modules; Cardiocap/5 and S/5 Light Monitors,
 - Ohmeda Trusat,
 - CARESCAPE™ modular monitors and S/5, AS/3 and CS/3 modular monitors with E-(P)RESTN or E-PSM(P) modules. CARESCAPE™ V100 with TruSignal® option, Responder 2000,
 - Handheld pulse oximeters (GE Ohmeda TuffSat®, 3775), tabletop pulse oximeters (3700 series, 3800, 3900, 3900P), respiratory gas monitors (4700, 5250 RGM), CARESCAPE modular monitors and S/5, AS/3, CS/3 modular monitors with M-OSAT module, CardiocapTM/5 with N-XOSAT option,
- Gwarancja 12-miesięcy,

